

Update on GS1 and UDI

23rd AHWP Annual Meeting Program


22 - 25 October 2018, Kuala Lumpur, Malaysia

Disclaimer



Neither GS1 nor its member organizations nor their staffs have *real* or apparent authority to speak for the regulatory authorities or grant exemptions. GS1 offers advisory services focused on the GS1 standard after a supplier's staff including its internal regulatory experts have determined the correct path to compliance. GS1 is a voluntary organization and its members have and must continue to determine their own course of action. GS1 provides recommendations. GS1 Global Office, GS1 member organizations and GS1 staff assume no liability for members actions taken upon its advice.







GS1 Standards





GS1 – a global non-for profit standard organisation





1 million

over 1 million companies worldwide use GS1 standards 150 countries

25 industries served across 150 countries

6 billion


Barcodes scanned more than 6 billion times per day globally **112 MOs**

112 Member Organisations around the world



A global system of standards to ensure visibility








GS1 standards for harmonized implementation of UDI





UDI system and GS1 system








UDI in AIDC terms







The HRI format shall follow the rules of the UDI Issuing Agency



UDI Example of label





H.E.L.P. Acetate Buffer pH 4.85

CA/GB Sodium acetate buffer solution for use ONLY with extracorporeal H.E.L.P. apheresis

Caution: Federal Law (U.S.) restricts this device to sale by or on order of a physician.

CA/FR: Solution tampon d'acetate de sodium destinée à une utilisation UNIQUEMENT avec aphérèse H.E.L.P. extracorporelle

sterile / stérile

Endotoxin-FREE and non-pyrogenic/ Ne contient pas d'endotoxines et non-pyrogène

SINGLE USE only, discard unused portion/ À USAGE UNIQUE seulement, jeter la portion inutilisée DO NOT add any additives/ NE PAS ajouter d'additifs

NOT for intravenous infusion/ NON adapté à une perfusion intraveineuse

ONLY USE if solution is clear and colourless/ UTILISER UNIQUEMENT si la solution est limpide et incoloré ONLY USE if container and connections are not damaged/ Ne pas utiliser si l'emballage et les

connections sont endommagées

Keep out of the reach of children/ Conserver la solution hors de portée des enfants

Sodium acetate x 3 H₂O 27.22 q/l Acetic acid 99% 6.82 q/l

DIN: 02373807



(01)0404696429016 (17)170228(11)140304(10)03502

Manufacturer:

B BRAUN

B. Braun Avitum AG 34209 Melsungen Germany

4 x 3000 ml









Read



Article no.: 4113

Batch no.: 0350214

Manuf. date: 2014-03-04



Expiry date: 2017-02-28

Canadian Distributor:

Chief Medical Supplies Ltd. 411-19th Street S.E. Calgary, Alberta T2E 6J7

US Distributor:

B. Braun Medical Inc. Bethlehem.PA 18018-3524

Production site:

B. Braun Avitum AG Kattenvenner Str..32 49219 Glandorf, Germany

Made in Germany

Production Identifier (PI)

"Dynamic" portion

Application Identifiers (e.g. serial, lot number & expiry date)




Device Identifier (DI)

"Static" portion

GTIN (product identifier)

NEW: EU specificity: "Basic UDI-DI" for illustration only








NEW GS1 Standard: "Global Model Number" to support implementation of the Basic UDI-DI



- GS1 Company Prefix + a model reference
- Model reference number: manufacturers internal numbers, can be alphanumeric



In healthcare:

- data title: "BUDI-DI"
- for medical device registration
- independent of packaging
- never used in a data carrier




 GS1 implementation guide to be drafted in 2018-2019, based on the EU Commission guidelines



UDI Databases – USA and EU








Managing data and global standards: Global Data Synchronisation Network (GDSN)





Manufacturers are able to provide data to all UDI databases and their customers (hospitals, distributors, wholesalers, GPOs) simultaneously, with a single connection

Mapping of GDSN data elements with EUDAMED will start when technical specifications will be released by the EU Commission



Time for change is NOW ...







Developments and requirements for UDI







GS1 key to implement requirements on medical devices



Accredited as UDI Issuing Agency by the US FDA

99% of medical devices identified with GTIN in Japan UDI assigning entities listed in the EU MDR



£3 million on average saved each year in every NHS hospital in England



91,8% of devices identified with GTIN in Turkey

Turkish National Drug and Medical Device Databank (TITUBB)

Mandated by ANMAT for traceability of certain devices in Argentina



The benefits of UDI





Manufacturers:

- cost optimization
- data synchronization
- process efficiency



Hospitals:

- adequately identified medical devices
- single and integrated system of information management



Regulators:

- higher levels of market surveillance
- more efficient adverse event reports and quicker recall also across borders



The need to align on a global UDI framework



 UDI is very beneficial. It is crucial that regulators around the world align on the global framework - IMDRF Guidelines and ensure consistency when setting-up regional or national UDI system

This will ensure:

- highest levels of patient safety beyond borders
- harmonized identification systems for medical devices globally



34th Global GS1 Healthcare Conference





30 October to 1 November 2018 - Bangkok, Thailand

- Traceability, UDI and global regulatory developments
- Success stories of implementations from manufacturers, wholesalers and hospitals – experiences & benefits
- Patient safety and quality of care how to improve those
- Awards for best hospital implementations
- And much more...

SAVE THE DATE!



http://www.gs1.org/healthcare/udi







Contact Details

Géraldine Lissalde-Bonnet
Director Public Policy
GS1 Global Office, Brussels
E g.lissalde@gs1.org
W www.gs1.org/healthcare